

OVERVIEW

PURPOSE: To develop a proper methodology to link flies to pathogen transmission via ingestion of infected feces

METHODS: Fly gut extract samples were analyzed via LC-MS/MS

RESULTS: Urobilin and urobilinogen were detected at a retention time of 6.5 min and 8.2 min respectively only in flies that fed on feces

INTRODUCTION

- No clear link stating flies transmit pathogens via ingestion of fecal matter.
- Most analysis are done via molecular biology techniques such as DNA sequencing or bacterial titers
 - Cumbersome
- No compounds found in literature that tie fly consumption of feces to pathogen transmission

METHODS

- Flies were fed based on control group
- Flies were killed, dissected, and DNA extractions were performed
- 50 μ L aliquot of organic layer from DNA extraction was evaporated under nitrogen
- Resuspend in 50 μ L of 1:1 methanol:water solution
- Vortex for 10 min
- Separate using Agilent 1100 HPLC system
- Identification on Thermo Fisher LTQ XL™ Linear Ion Trap Quadrupole




Figure 1: LC-MS instrument set up




Figure 2: Structures for urobilin and urobilinogen




Figure 3A: Extracted ion chromatogram and mass spectra for urobilin




Figure 3B: Extracted ion chromatogram and mass spectra for urobilinogen




Figure 4: Extracted ion MS/MS chromatograms (XIC) and MS/MS spectra for fly samples (absolute intensities).




Figure 5: Scatterplot comparison of LC MS/MS 6.3 – 6.6 and 8.2 min peak area data for all tissue and fecal controls, as well as all experimental and wild flies

from Feces in Fly Guts

Christine Picard and Nick Manicke

Indiana University-Purdue University Indianapolis
Indiana University-Purdue University Indianapolis

ULTS

CONCLUSIONS




Figure 6A: Xbar control chart for the urobilin positive control samples (RT ~6.5 min) showing overall process variation




Figure 6C: Xbar control chart for the urobilinogen positive control samples (RT ~8.2 min) showing overall process variation




Figure 6B: Range control chart for the urobilin positive control samples (RT ~6.5 min) showing overall process variation




Figure 6D: Range control chart for the urobilinogen positive control samples (RT ~8.2 min) showing overall process variation





Figure 7: Individual moving range control charts of the peak areas for the urobilinoid standard showing overall instrument variation

Conclusions:

- LC MS/MS is a good qualitative test to detect fecal urobilinoids
- No false positives or negatives have been observed
- According to the Kruskal-Wallis test, there was a significant effects of treatments at 6 min (chi-squared = 72.8, df = 3, P < 0.0001) and 8 min (chi-squared = 78.3, df = 3, P < 0.0001)

Potential Problems:

- Many factors can impact the urobilinoid signal intensity
 - Feces consumed by the fly
 - Exposure of samples to light and air
 - Gender of the fly
- Overall instrument variation
- Difficult to get "pure" urobilin standards
- Difficult to fully resolve urobilinoid isomers
 - Second peak at ~8.2 min helped with identification of these
- High variation among positive "control" samples
 - Standard was used to help with this

Future Directions:

- Apply method using multiple species of filth flies and other coprophagous insects
- Combine the method with microbial culturing and sequencing methods
- Combine the method with vertebrae sequencing methods to identify the source of pathogens

REFERENCES

[1] CG Owings, C Skaggs, W Sheriff, N Manicke, CJ Picard. Chemical assay for the detection of vertebrate fecal metabolites in adult blow flies (Diptera: Calliphoridae). Environmental Entomology

ACKNOWLEDGEMENTS

We would like to thank the Indiana University- Purdue University Indianapolis Department of Chemistry and Biology for allowing us to perform these experiments